

Technical data

MCC 250K100M

Casing material: Steel lacquered

Dimensions of display device (W×D×H): 200 mm x 130 mm x 60 mm

Dimensions of weighing plate (W×D×H): 502 mm x 380 mm x 410 mm

Level indicator: yes

Material weighing plate: Plastic

Overall dimensions mounted (W×D×H): 625 mm x 990 mm x 985 mm

Weighing surface (W×D): 390 mm x 360 mm

Adjustment at the location of installation: 961-249

DAkKS Certificate: 963-129

Verification: 965-129

Display digit height (large): 25 mm

Technical data about verification

Minimum load [Min]: 2 kg

Periodic Verification for Italy: 969-112

Verification approval conforming to 2014/31/EU: yes

Verification class: III

Verification value [e]: 100 g

Battery: 6×1.5 V AA

Charging time: 12 h

Input voltage: 100 V - 240 V 50 / 60 Hz

Mains adapter: Mains adapter external

Mains adapter/adaptor included: EURO CH UK

Operating time: 20 h

Operating time (Backlight off): 40 h

Operating time (Backlight on): 20 h

BMI function: yes

Data hold function: yes

Vibration-free weighing (Animal weighing program): yes

Readability [d]: 100 g

Recommended adjusting weight: 200 kg (M1)

Repeatability: 100 g

Resolution: 2500

Stabilisation time under laboratory conditions: 3 s

Tare range: 249,900 kg

Warm up time: 10 min