



GUARANTEED  
**3**  
YEARS

## ASPIRATORE PROFESSIONALE PROFESSIONAL ASPIRATOR

### DATI TECNICI VERSIONE POMPA SINGOLA

BASSO FLUSSO .....	Alto vuoto
PORTATA D'ARIA .....	15 lt/min
LIVELLO VUOTO REGOLABILE .....	0 ÷ - 0,85 bar (0 ÷ - 85 kPa) reg.
CICLO DI UTILIZZO .....	Continuo
RUMOROSITÀ (a 1 mt) .....	55 dBA
PESO .....	3,5 kg

### DATI TECNICI VERSIONE DOPPIA POMPA

BASSO FLUSSO .....	Alto vuoto
PORTATA D'ARIA .....	22 lt/min
LIVELLO VUOTO REGOLABILE .....	0 ÷ - 0,85 bar (0 ÷ - 85 kPa) reg.
CICLO DI UTILIZZO .....	Continuo
RUMOROSITÀ (a 1 mt) .....	65 dBA
PESO .....	4,5 kg

### ACCESSORI IN DOTAZIONE

VASO DI RACCOLTA DA 1000cc IN POLICARBONATO, AUTOCLAVABILE, con valvola di troppo pieno  
SACCA MONOUSO DA 1000cc  
FILTRO ANTIBATTERICO 99% IDROFOBICO + 1 RICAMBIO  
CANNULA STERILE MONOUSO  
REGOLATORE DI FLUSSO MANUALE, STERILE, MONOUSO  
SET TUBI IN SILICONE ATOSSICO STERILIZZABILE  
CAVO DI ALIMENTAZIONE

### TECHNICAL DATA SINGLE PUMP VERSION

LOW FLOW .....	High Vacuum
AIR FLOW .....	15 lt/min
ADJUSTABLE VACUUM LEVEL .....	0 ÷ - 0.85 bar (0 ÷ - 85 kPa) adjus.
OPERATING CYCLE .....	CONTINUOUS
NOISE LEVEL (at 1 mt) .....	55 dBA
WEIGHT .....	3.5 kg

### TECHNICAL DATA DOUBLE PUMP VERSION

LOW FLOW .....	High Vacuum
AIR FLOW .....	22 lt/min
ADJUSTABLE VACUUM LEVEL .....	0 ÷ - 0.85 bar (0 ÷ - 85 kPa) adjus.
OPERATING CYCLE .....	CONTINUOUS
NOISE LEVEL (at 1 mt) .....	65 dBA
WEIGHT .....	4.5 kg

### EQUIPPED WITH FOLLOWING ACCESSORIES

AUTOCLAVABLE POLYCARBONATE JAR 1000cc with safety valve (overflow protection)  
DISPOSABLE SUCTION LINER 1000cc  
99% ANTIBACTERIAL HYDROPHOBIC FILTER + 1 SPARE  
STERILE, DISPOSABLE CANNULA  
STERILE, DISPOSABLE, MANUAL FLOW REGULATOR  
SET OF ATOXIC STERILIZABLE SILICONE TUBES  
POWER CORD

CE0434

(dir. 93/42/EEC)

PCS.

2

MASTER CARTON



DISPOSITIVO MEDICO PER ASPIRAZIONE  
CHIRURGICA PROFESSIONALE

MEDICAL DEVICE FOR PROFESSIONAL SURGICAL  
ASPIRATION

