

LHH | LHH LED
LHH024803



LHH 10 W105 WH 50W GU5.3 GB

LHH is a medical examination luminaire with excellent light output, mounted on a flexible self-balancing arm with frictionless horizontal and vertical movement.



Electrical data

Voltage from (V)	230
Maximum voltage (V)	240
Maximum frequency (Hz)	60
Minimum frequency (Hz)	50
Total consumption (W)	50

Technical data

Maximum ambient temperature (°C)	25
IP classification	20

Optic

Optic	SYM
Luxlevel at 1m	4330

Light source

Number of lamps	1
Lamp power (W)	50
Lamp base	GU5.3
Light source	QR-CBC
Colour rendering and colour temperature	933

Body

Colour of body	White (RAL9016)
Body material	Steel

Termination

Plug type	UK
Cable length (m)	2,9

Dimensions

Arm length (cm)	105
-----------------	-----

LHH | LHH LED
LHH024803



LHH 10 W105 WH 50W GU5.3 GB

LHH is a medical examination luminaire with excellent light output, mounted on a flexible self-balancing arm with frictionless horizontal and vertical movement.



Drawings

