



DispozABLE Spacer

Technical specifications

Tradition and Innovation - Visionary thinking and a fascination with technology have guided us throughout our history to develop innovative products of outstanding reliability. Anticipating trends to improve the quality of life.

DispozABLE Spacer

Technical specifications

DispozABLE™ Spacer

DispozABLE Spacer

Dimensions	ø84mm x 107mm
Weight	15g
Volume	227ml
Resistance to flow	Free flow

Performance Data

Fine particle dose (salbutamol):	
Delay 0 secs =	54.5mcg (SD 4.8)
Delay 1 secs =	37.7mcg (SD 7.6)

Product	Part Code
DispozABLE Spacer (box of 10)	3609550

Product information

DispozABLE Spacer has been designed to meet a specific market need. It is a low cost, fast response spacer for delivering medication promptly. It is designed to be used straight from the box without pre-washing in an emergency situation. For example in an ambulance or school. DispozABLE Spacer is to be used as a single use spacer. This simple intuitive design allows non-professionals to assemble it quickly. Clear instructions are included on the spacer, adding to the element of ease.

- Simple to use
- Effective
- For use in an emergency situation: ambulance or school
- Easy to assemble
- No mask required
- Can be used with tidal breathing technique, assisted by Flo-Tone.