

BM3 AND BM5 PORTABLE MULTIPARAMETER VET MONITORS



33779

- 33719 BM3 VET MONITOR - 7" TFT display
 - 33779 BM5 VET MONITOR - 10.4" TFT display
- Versatile, portable, lightweight, BM monitors are designed to be used in veterinary environment. Compact design and handle make it easy to move. Colour monitor displays all vital signs in different colours for easy recognition. A built in printer offers hard copy documentation of all measured values. Alarms and settings are easily adjustable. See all specifications and accessories at page 245-246.

- 33741 EXACON RECTAL TEMPERATURE SENSOR - 2.75 m cable

ECG+SpO₂+NIBP
+RESP+TEMPERATURE
+2IBP (only 33779)



3-lead ECG clip



SpO₂ sensor

HANDHELD VET PULSE OXIMETER



- 34343 OXY 100 VET PULSEOXIMETER
- OXY 100 Handheld PulseOximeter can measure oxygen saturation (SpO₂) and pulse Rate (PR) by placing the sensor on the certain part of the animal, such as tongue, ear, leg or tail. It also provides plethysmogram, bar-graph and Perfusion Index display.
- Lightweight, small and easy to handle
 - Data storage with trend review
 - Two options for measurement on: thin tissue: such as tongue, ear, thick tissue: such as leg, tail.
 - Up to 384 hours of data memory (with 1/2/4/8 sec. recording interval)
 - Audible & Visual alert for over-limits (limits are adjustable).
- Oximeter is supplied with universal "Y" type SpO₂ probe. See technical specifications page 251.
- 34346 UNIVERSAL "Y" type SpO₂ probe - spare includes 3 types of probe adapters (big clip, small clip, rubber wrapper).



34346

OXY 100 is equipped with a special probe for animals monitoring. The universal "Y" type SpO₂ probe can be used in combination with the provided three types of adapters: big clip, small clip and rubber wrapper. Each of these adapters allows a different measuring respectively on tongue, ear, leg or tail.

PREVITNT COHESIVE NON WOVEN ELASTIC BANDAGES

PREVITNT NON WOVEN ELASTIC BANDAGES

For fixing dressings. Ideal for supporting joints, treating light fractures. This bandage does not require clips because it remains well-attached to itself. Particularly suitable for use on joints, for sports and veterinary use. Composition: 85% polypropylene - 15% spandex and latex. Extensibility: approx 100%.



COHESIVE COLOUR NON WOVEN BANDAGE				
Size	Blue	Red	5 Colours*	Minimum order
4.5 m x 7.5 cm	34734	34735	34736	box of 10
4.5 m x 10 cm	34738	34739	34740	box of 10

*Mix of 2 bandages x 5 colours (blue, red, green, yellow, skin)