

Technical Specifications

1. Model	AL-01
2. Device Class	Class IIa Medical Device
(Medical Device Directive No. 93/42/EEC)	
3. EN ISO 10079-1	High vacuum/Low air flow
4. Power Supply	AC220V \pm 10% - 50Hz \pm 2%
5. Input Power	180VA
6. Fuse	F 1 x 1.6AL 250 V
7. Maximum vacuum Pressure limit (without jar)	-0.75bar -75kPa -563mmHg
8. Noise	\leq 55dB (A)
9. Liquid collection jar	1000 ml
10. Weight	2.41 kg
11. Dimensions	240x190x130 mm
12. Operation system	Continuous operation

- The device must not be used in environments with flammable and explosive gas.