

# VacuAide 7310 Compact Suction Unit



Discreet  
suction on  
the go

The VacuAide 7310 is one of the smallest and strongest suction units on the market. With its lightweight and small size and featuring a carry bag and an enhanced battery life; this suction unit is designed to fit the active lifestyles of today.

## VacuAide 7310 - Compact Suction Unit



### Features

- » Vacuum adjustment allows for 50 – 550 mm Hg and a free flow of 27 litres/minute
- » Enhanced battery of life for up to 1 hour of operation
- » Meets ISO1007901:1999 high vacuum, high flow classification
- » Easy to use LED bar gauge which provides precise measurements of vacuum – especially important for paediatric patients
- » Compact and lightweight design; weighing only 1.51kg



### Specifications

Code  
7310PR-D

Length  
7.25"

Width  
6.75"

Height  
7.25"

Weight  
3.37lb

Vacuum Range  
50 – 550 mmhg

Air Flow  
27 l/min

Electrical Requirements  
100 – 240 VAC, 50/60 Hz

Internal Battery Run Time  
60 minutes

Warranty  
2 years unit, 6 months battery

### Bottle design

- » Reusable autoclavable 725ml bottle
- » In-built unique float valve bacteria filter in the suction bottle prevents backflow into the unit



### Portable suctioning

- » Supplied with its own carry bag designed for discreet transportation and suctioning
- » Travel-friendly design meets RTCA/DO-160D Aircraft Standard for Environmental Conditions and Test Procedures for Airborne Equipment

