

Able Spacer Technical specifications

Tradition and Innovation - Visionary thinking and a fascination with technology have guided us throughout our history to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.



Able Spacer Technical specifications



Able Spacer	
Dimensions	155 x 55 x 66mm
Weight	62g
Volume	210ml
Resistance to flow	10mm of water at 30 litres/minute

Product	Part Codes
Able Spacer in a box	3609120
Able Spacer with small whistle mask in box	3609121
Able Spacer with medium whistle mask in box	3609122
Able Spacer with large mask in box	3609123
Able Spacer Masks	
Able Spacer Mask Small Infant (pack of 1)	3609150
Able Spacer Mask Medium Child (pack of 1)	3609151
Able Spacer Mask Large Adult (pack of 1)	3609152
Able / A2A Small Mask (box 50)	3607125
Able / A2A Medium mask (box 50)	3607126

Product information

Able Spacer is a popular valved inhalation chamber that facilitates the use of pMDI medication. Able Spacer and Mask are made from a low-static transparent plastic, which incorporates anti-microbial silver ion technology; this protects against contamination and inhibits microbial growth. This improved version of Able Spacer, manufactured in the UK, now features a two-way valve; operating at less than 10L/min inspiratory flow. Recommended for children and COPD patients.

- Inhibits microbial growth
- Low-static properties
- Two-way valve
- Visible valve movement
- 210ml, consistent with optimal performance
- Supplied in an anti-static bag
- Includes a tethered mouthpiece cap
- Available with or without mask

Able Spacer is manufactured in accordance with all appropriate EU regulations.

© CLEMENT CLARKE INTERNATIONAL Ltd,